

# MasterCare

Patient Equipment, Inc.

*"ERGONOMIC AND ECONOMIC PATIENT CARE SOLUTIONS"*

---

## Recumbent Bath

### Architectural Specifications & Rough-In Information

#### Dimensions:

<b>Bath:</b>	Width	35 ¼" Inches
	Length	84" Inches
	Rail Height	34 ½" Inches
	Deck Height	38 ¾" Inches

#### Specifications:

<b>Back Flow Prevention:</b>	Air gap for fill spout. Note: The factory provided backflow prevention is the most widely accepted. Local codes may require changes, which is the responsibility of the local plumbing contractor.	
<b>BathAire<sup>(SM)</sup> Outlets:</b>	4 with 6 ports each for a total 24 warm air massage outlets.	
<b>Electrical:</b>	Provide separate 115 V-AC, 60 HZ, 20 Amp. circuit and receptacle with in-line ground fault interrupter (not supplied) in accessible location.	
<b>Installation Time:</b>	1 hour	<b>Inspection and Adjustments:</b> ½ hour
<b>Motor Specs:</b>	BathAire <sup>(SM)</sup> Motor: 115 V., 7.0 Amps, 1.3 H.P. Adjustable speed for maximum comfort and control.	
<b>Room Size:</b>	7' wide x 8' 10" long minimum recommended. Can fit in smaller room; call MasterCare with dimensions.	
<b>Switches:</b>	Membrane.	
<b>Thermometer:</b>	Analog temperature measurements.	
<b>Tub Fill/Shower:</b>	30 GPM thermoscopic mixing valve with separate fill and shower flow control valves.	
<b>Tub Shell:</b>	Sanitary gel-coated FRP (Fiberglass Reinforced Polymer).	
<b>Tub Water Use Gallonage:</b>	65 gallons	
<b>Waste:</b>	2" PVC schedule 40 pipe.	
<b>Water Supply:</b>	¾" Copper hot and cold. Minimum of 30 PSI for hot and cold required. Maximum temperature of hot water supply 115° F or code, whichever is less.	
<b>Weight:</b>	400 lbs.	

Recumbent Specs 10/25/06

---

2071 14<sup>th</sup> Avenue · PO Box 1435 · Columbus, NE 68601

Phone: 800.798.5867 · Fax: 402.563.9102

E-mail: [mastercarepeinc@frontiernet.net](mailto:mastercarepeinc@frontiernet.net)

Website: [www.mastercarebath.net](http://www.mastercarebath.net)